

## DNA extraction from linseed containing samples using SureFood® PREP Basic

February 2025

### 1. Sample preparation

Please prepare the samples in duplicates

### 2. Additional required equipment and materials

- micro balance for weighing the sample
- 2.0 ml reaction tubes (not provided with the kit)
- heating block (up to 65°C)
- micro centrifuge (up to 12,000 rpm)

### 3. Procedure

- Transfer 50 mg sample material into a 2.0 ml reaction tube
- Add 800 µl Lysis Buffer and 40 µl Proteinase K
- Mix briefly on a Vortex mixer
- Incubate on a heating block under continuously shaking for 30 min at 65°C, according to the User Manual step 2 Lysis of the basic material
- Centrifuge the sample lysate 1 min at 12,000rpm
- Place a Spin Filter (Code F, for pre-filtration) into a 2.0 ml Receiver Tube (Code R)
- Carefully transfer the supernatant directly onto the Spin Filter, without the oil phase or precipitated proteins
- Centrifuge the Spin Filter with the Receiver Tube 1 min at 12,000rpm
- After centrifugation discard the Spin Filter
- Add 400 µl Binding Buffer to the filtrate and mix
- Place a Spin Filter (Code S, for binding) into a new 2.0 ml Receiver Tube (Code R)
- The binding of the nucleic acids on a Spin Filter takes place in two steps
  1. Transfer 650 µl of the filtrate-binding buffer mixture to the Spin Filter
  2. Incubate at room temperature for 1 min
  3. Centrifuge the Spin Filter with the Receiver Tube 1 min at 12,000 rpm
  4. Discard the filtrate and place the Spin Filter back into the Receiver Tube
  5. repeat the steps 1. to 4.
- Continue with step 5 Purification of the bound nucleic acids from SureFood® PREP Basic