**CONGEN**

# **SureFast<sup>®</sup> Legionella 3plex**

Art. No. F5505  
100 rxn

## **User Manual**

**July 2019**

 **Inhalt**

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 **Content**

1	General Information .....	7
1.1	Description .....	7
1.2	Limit of Detection .....	7
1.3	DNA-preparation .....	7
1.4	Kit components and storage .....	7
1.5	Additionally required equipment and materials .....	7
1.6	Setup .....	8
1.7	Detection channel Set-up .....	8
2	Qualitative Analysis .....	9
2.1	Protocol .....	9
2.1.1	Preparation of the master-mix .....	9
2.1.2	Preparation of the real-time PCR-mix .....	9
2.2	Interpretation of results .....	10
3	Further Information .....	10
3.1	Product Information .....	10
3.2	Technical Support .....	10

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## 1 General Information

### 1.1 Description

The SureFast® Legionella 3plex is a real-time PCR for the direct, qualitative detection and differentiation of specific DNA sequences of *Legionella* spp. and *Legionella pneumophila*.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX and Cy5 at the same time. The technical validation of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Bio-Rad CFX96, Roche LightCycler® 480 II, Roche cobas z 480 Analyzer, Applied Biosystems 7500, Qiagen Rotor-Gene Q and R-Biopharm RIDA®CYCLER.

### 1.2 Limit of Detection

The SureFast® Legionella 3plex real-time PCR has a limit of detection of  $\leq 5$  DNA copies. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Aqua is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP Aqua Art. No. F1023)
- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

**1.6 Setup**

	<b>Blockcycler &amp; R-Biopharm RIDA®CYCLER</b>	<b>Rotorcycler &amp; LightCycler® 480 II</b>
Initial Denaturation (HOLD) Cycles	1 min, 95°C 45	1 min, 95°C 45
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 60°C	10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.7 Detection channel Set-up**

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent Mx3005P</b>	<i>Legionella</i> spp.	FAM	+	
	IAC	HEX	+	
	<i>Legionella pneumophila</i>	Cy5	+	
<b>Applied Biosystems 7500</b>	<i>Legionella</i> spp.	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
	<i>Legionella pneumophila</i>	Cy5	None	
<b>Bio-Rad CFX96</b>	<i>Legionella</i> spp.	FAM	+	
	IAC	VIC	+	
	<i>Legionella pneumophila</i>	Cy5	+	
<b>R-Biopharm RIDA®CYCLER</b>	<i>Legionella</i> spp.	green	+	
	IAC	yellow	+	
	<i>Legionella pneumophila</i>	red	+	
<b>Qiagen Rotor-Gene Q</b>	<i>Legionella</i> spp.	green	+	
	IAC	yellow	+	
	<i>Legionella pneumophila</i>	red	+	
<b>Roche LightCycler® 480 II</b>	<i>Legionella</i> spp.	465-510	+	
	IAC	533-580	+	
	<i>Legionella pneumophila</i>	618-660	+	
<b>Roche cobas® z 480 Analyzer</b>	<i>Legionella</i> spp.	465-510	+	
	IAC	540-580	+	
	<i>Legionella pneumophila</i>	610-670	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. For the analysis of enrichments additional controls are recommended: zero control (sample before enrichment) and medium control. The reaction mix contains an internal amplification control (IAC) per reaction.

#### Reactions needed for the qualitative *Legionella* spp. and *Legionella pneumophila* detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

#### Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

**2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

*Legionella* spp. DNA is detected in the FAM-channel and *Legionella pneumophila* DNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-value  $\leq 2$  compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific *Legionella* spp. or *Legionella pneumophila* PCR assay.

Result in the respective channel			Interpretation
FAM channel <i>Legionella</i> spp.	Cy5 channel <i>Legionella pneumophila</i>	VIC/HEX channel IAC	
<b>positive</b>	negative	<b>positive</b>	<i>Legionella</i> spp. DNA detected
<b>positive</b>	<b>positive</b>	<b>positive</b>	<i>Legionella pneumophila</i> DNA detected
negative	negative	negative	invalid

**3 Further Information**

**3.1 Product Information**

- Detailed information about setup of several real-time PCR devices (Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request

**3.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).